
































































SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
1.		CE marking	Indicates that the product meets the requirements for accreditation and market surveillance relating to the marketing of the products.	Regulation (EU) 2017/745 or MDD 93/32/EEC
2.		CE marking identification number by TUV SUD Product Surveillance GmbH Germany	Indicates that the product meets the requirements for accreditation and market surveillance relating to the marketing of the products.	Regulation (EU) 2017/745 or MDD 93/32/EEC
3.		UK CA marking	UK product marking that is used for goods being placed on the market in Great Britain (England, Wales and Scotland).	UK MDR 2002
4.		Authorized Representative in UK	Indicates the authorized representative in UK.	UK MDR 2002
5.		Manufacturer	Indicates the medical device manufacturer.	ISO 15223-1 Reference no. 5.1.1
6.		Authorized representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1 Reference no. 5.1.2
7.		Date of manufacture	Indicates the date when the medical device was manufactured.	ISO 15223-1 Reference no. 5.1.3
8.		Use by date	Indicates the date after which the medical device is not to be used.	ISO 15223-1 Reference no. 5.1.4
9.		Batch code	Indicates the manufacturers batch code number or lot number for identification purpose.	ISO 15223-1 Reference no. 5.1.5
10.		Catalogue number	Indicates the medical device catalogue reference number for identification purpose.	ISO 15223-1 Reference no. 5.1.6
11.		Serial number	Indicates the medical device serial number for identification purpose.	ISO 15223-1 Reference no. 5.1.7
12.		Sterile	Indicates that the product has been subjected to a sterilization process.	ISO 15223-1 Reference no. 5.2.1

SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
13.		Sterilized using ethylene oxide	Indicates that the product has been sterilized using ethylene oxide.	ISO 15223-1 Reference no. 5.2.3
14.		Sterilized using irradiation	Indicates that the product has been sterilized using irradiation.	ISO 15223-1 Reference no. 5.2.4
15.		Non-sterile	Indicates a device that has not been subjected to a sterilization process.	ISO 15223-1 Reference no. 5.2.7
16.		Do Not Resterilize	Indicates a medical device that is not to be resterilized.	Curas created symbol
17.		Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damage or opened.	ISO 15223-1 Reference no. 5.2.8
18.		Keep away from sunlight	Indicates that the product needs protection from light sources.	ISO 15223-1 Reference no. 5.3.2
19.		Keep dry	Indicate that the product needs to be protected from moisture / water.	ISO 15223-1 Reference no. 5.3.4
20.		Lower limit of temperature	Indicates the lower limit of temperature to which the medical device can be safely exposed.	ISO 15223-1 Reference no. 5.3.5
21.		Upper limit temperature	Indicates the upper limit of temperature to which the medical device can be safely exposed.	ISO 15223-1 Reference no. 5.3.6
22.		Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	ISO 15223-1 Reference no. 5.3.7












SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
23.		Humidity limitation	Indicates the range of humidity to which the medical device can be safely exposed.	ISO 15223-1 Reference no. 5.3.8
24.		Atmospheric pressure limitation	Indicates the range of atmospheric pressure to which the medical device can be safely exposed.	ISO 15223-1 Reference no. 5.3.8
25.		Biological risks	Indicates that there are potential biological risks associated with the medical device.	ISO 15223-1 Reference no. 5.4.1
26.		Do not re-use	Indicates a medical device that is intended for one use, or on a single patient during a single procedure.	ISO 15223-1 Reference no. 5.4.2
27.		Consult instructions for use	Indicates the need for the user to consult the instructions for use.	ISO 15223-1 Reference no. 5.4.3
28.		Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot for a variety of reasons, be presented on the product itself.	ISO 15223-1 Reference no. 5.4.4
29.		Contains or presence of natural rubber latex	Indicates the presence of natural rubber or dry natural rubber latex as a material of construction within the medical device or the packaging of a medical device.	ISO 15223-1 Reference no. 5.4.5
30.		Latex-free	Indicates that natural rubber latex was not used in the manufacturing of the product, its container, or its packaging.	ISO 15223-1 Reference no. 5.4.5 Reference Annex B for the general prohibition symbol and negation symbol
31.		Non-pyrogenic	Indicates that the medical device is non-pyrogenic.	ISO 15223-1 Reference no. 5.6.3
32.		Phthalate (DEHP)-containing device	Contains or presence of phthalates bis (2- ethylhexyl) phthalates (DEHP).	BS EN 15986 EN 15986 (E) clause 4.2 and Annex A

SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
33.		RoHS Complaint	Indicates that this product is in compliance to the European Union's Directive 2002/95/EC, the Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment.	European Union's Directive 2002/95/EC
34.		Medical Device	Indicates that the device is a medical device as defined in MDR 2017/745.	Regulation (EU) 2017/745
35.		Phthalate (DEHP)-free	Indicates that the product do not contain phthalate (DEHP).	Curas created symbol
36.	<b>EN 455</b> (Parts 1, 2 & 3)	European standard for medical gloves	European standard for medical gloves.	Curas created symbol
37.	<b>AQL 1.5</b>	AQL of 1.5 for pinholes	Acceptable quality levels of 1.5 for pinholes.	Curas created symbol
38.		Gloves giving protection from chemicals and micro-organisms	European standard: Protective gloves against chemicals and microorganisms. Microorganism resistant.	Standard EN374: 2016
39.		Protection against mechanical hazards	European standard: Protective gloves against mechanical risks.	Standard EN388: 2016
40.		Bacterial Filtration Efficiency (BFE)	Indicates that the product conform to BFE Type IIR test.	Standard EN 14683:2019
41.		PTFE Filter (3 Plus 1™)	Indicates that the product is filtered by 3 layers + 1 layer coating of polytetrafluoroethylene material.	Curas created symbol
42.		Fluid resistance test PAS 29:1999	Indicates fluid resistance test for a minimum of 4 hours (water at 35°C±3°C) has been conducted.	Curas created symbol











SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
43.		Personal Protective Equipment, Category 1	Personal Protective Equipment, European Directive 89/686/EEC, Category I: Protects from minor risks only.	European Directive 89/686/EEC
44.		Polyvinyl chloride-free	Indicates that the product do not contain polyvinyl chloride.	Curas created symbol
45.		Do not flush	Indicates that this product should not be flushed down the toilet bowl.	Curas created symbol
46.		Period-after-opening	Indicates that product should be used within 6 months after the seal is opened.	Curas created symbol
47.		Powder	Indicates that this product contains traces of powder.	Curas created symbol
48.		Powder-free	Indicates that this product does not contain any traces of powder.	Curas created symbol
49.		Macerator friendly	Indicates that this product can be disposed through maceration.	Curas created symbol
50.		Macerator friendly	Indicates that pulp component inside this product can be disposed through maceration.	Curas created symbol
51.		Macerator friendly	Indicates that this product can be disposed through maceration. Performance is subject to individual macerator. Seek for technical advice when necessary.	Curas created symbol
52.		Pyrogenic-free	Indicates that this product does not contain any pyrogenic material.	Curas created symbol
53.		Male connector	Indicates that the connector can fit all common catheters.	Curas created symbol
54.		For inflation of balloon catheter only	Indicates that this product is intended for catheter balloon inflation only.	Curas created symbol











SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
55.		Not for injection	Indicates that this product is not for intramuscular injection.	Curas created symbol
56.		Safe for food contact	Indicates that the material used in the product is considered safe for food contact.	Curas created symbol
57.		Not for oily food	Indicates that the material used in the product is not for oily food.	Curas created symbol
58.		Non-return valve-free	Indicates that this product does not come with non-return valve.	Curas created symbol
59.		Non-return valve	Indicates that this product comes with non-return valve.	Curas created symbol
60.		Kink-resistant	Indicates that this product is resistant to kinking and will return to its normal shape.	Curas created symbol
61.		Needle free sampling port connector	Indicates that this product comes with needle-free sampling port connector for sampling procedure.	Curas created symbol
62.		Soap resistant up to 2 hours	Indicates that the product is resistant to soap for maximum 2 hours.	Curas created symbol
63.		Polymer coated	Indicates that the product is coated with polymer.	Curas created symbol
64.		Fully textured	Indicates that this product is with fully textured surface.	Curas created symbol
65.		Stainless steel 304	Indicates that this product is produced with stainless steel grade 304.	Curas created symbol












SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
66.		Autoclavable	Indicates that this product can withstand the high autoclave temperature.	Curas created symbol
67.		Ward Pack	Indicates that the product are packed in inner bags to protect from dust and contamination. It allows easy dispensing suitable for ward distribution.	Curas created symbol
68.		Clean Pack	Indicates that the product is packed individually in PE material.	Curas created symbol
69.		Room Pack	Indicates that the product is packed in inner bag that is suitable for individual room distribution.	Curas created symbol
70.		Box Pack	Indicates that the product comes in a box that allows easy dispensing.	Curas created symbol
71.		Flat Pack	Indicates that the product are folded in flat based packaging that is suitable for easy dispensing.	Curas created symbol
72.		Roll Pack	Indicates that the product is packed in roll form that is easy to pull off.	Curas created symbol
73.		Perfume-free	Indicates that this product does not contain any perfume.	Curas created symbol
74.		Colour free	Indicates that this product does not contain any colouring.	Curas created symbol
75.		Low lint	Indicates that this product is produced with low lint material.	Curas created symbol
76.		Lint-free	Indicates that this product does not contain any lint material.	Curas created symbol












SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
77.		High-density polyethylene	Indicates that the material/plastic used in this product is high-density polyethylene (HDPE).	Industrial Standard
78.		Low-density polyethylene	Indicates that the material/plastic used in this product is low-density polyethylene (LDPE).	Industrial Standard
79.		Polypropylene	Indicates that the material/plastic used in this product is polypropylene.	Industrial Standard
80.		Others	Indicates that the product is made of any combination of 1-6 type of plastic or another less commonly used plastic.	Industrial Standard
81.		Cardboard	Indicates that the packaging material is made of corrugated fiberboard (cardboard).	Industrial Standard
82.		Other paper	Indicates that the packaging material is made of non-corrugated fiberboard (paperboard).	Industrial Standard
83.		LDPE + 10% EVA	Indicates that this product is made of combination of LDPE and 10% of ethylene vinyl acetate.	Curas created symbol
84.		Slit cut	Indicates that this product comes with slit cut for easy opening.	Curas created symbol
85.		Flexi-Funnel™	Indicates that this product comes with Flexi-Funnel™.	Curas created symbol
86.		Medium Flexi-Funnel™	Indicates that this product comes with medium (M) size Flexi-Funnel™.	Curas created symbol
87.		Chlorine-free	Indicates that this product does not contain any chlorine.	Curas created symbol























SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
88.		Antimicrobial	Indicates that this product contain antimicrobial agents.	Curas created symbol
89.		Hypoallergenic adhesive	Indicates that this product comes with hypoallergenic adhesive.	Curas created symbol
90.		Compostable	Indicates that this product has been certified industrially compostable according to the European standard EN 13432.	Curas created symbol
91.		Biodegradable	Indicates that this product is capable of being decomposed by bacteria or other living organisms and thereby reduce pollution.	Curas created symbol
92.		Natural fibre	Indicates that this product is made from natural fibre material.	Curas created symbol
93.		Eco-friendly	Indicates that this product is not harmful to the environment.	Curas created symbol
94.		Recycled Packaging Material	Indicates that the packaging material is recycled.	Curas created symbol
95.		Recycled Packaging Material (percentage context)	Indicates that the packaging material is recycled with percentage context.	Curas created symbol
96.		Do not eat	Indicates that this product should not be eaten.	Curas created symbol
97.		Keep out of reach from children	Indicates that this product should be kept away from children.	Curas created symbol

SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
98.		Do not throw into basin	Indicates that this product should not be thrown into the basin.	Curas created symbol
99.		Super Absorbent	Indicates that the product absorbs water.	Curas created symbol
100.		Cross Hanger	Indicates that the product can be hanged on the bottom outlet of drainable urine bag.	Curas created symbol
101.		Bulb-pump	Indicates that the product has a bulb pump for cleaning/irrigation.	Curas created symbol
102.		SAP by Curas™	Indicates that Super Absorbent Polymer (SAP) is included in the product.	Curas created symbol
103.		Hydrophobic	Indicates that the product is water-repellent.	Curas created symbol
104.		Ultrasonic Seal	Indicates that the product is tightly sealed by ultrasonic heat seal.	Curas created symbol
105.		Light duty	Indicates that the product is meant for light duty loading.	Curas created symbol
106.		Medium duty	Indicates that the product is meant for medium duty loading.	Curas created symbol
107.		Heavy duty	Indicates that the product is meant for heavy duty loading.	Curas created symbol

SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
108.		Z folding	Indicates that the product is Z folded for easy dispensing.	Curas created symbol
109.		Z+1/2 folded for easy dispensing	Indicates that this product is Z and 1/2 folded for easy dispensing.	Curas created symbol
110.		Safe for use in the freezer	Indicates that the product is safe to be stowed in the freezer.	Curas created symbol
111.		Star-shaped kink-resistant tube	Indicates that this product comes with star-shaped kink-resistant tube.	Curas created symbol
112.		Safe for microwave	Indicates that the product is safe to be used in a microwave oven.	Curas created symbol
113.		Safe for oven	Indicates that the product is safe to be used in an oven.	Curas created symbol
114.		Wash at 30°C	Indicates that this product can be washed in water temperature up to 30 degrees Celsius.	Curas created symbol
115.		Half folded for easy dispensing	Indicates that this product is half folded for easy dispensing.	Curas created symbol
116.		Formaldehyde-free	Indicates that this product does not contain any formaldehyde.	Curas created symbol
117.		Open	Indicates the direction to open the product.	Curas created symbol
118.		High Capacity Cross Valve	Indicates that this product comes with High Capacity Cross Valve for efficient emptying.	Curas created symbol

SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
119.		Copolymer polyethylene	Indicates that this product contains copolymer polyethylene.	Curas created symbol
120.		Free Flow Cross Valve	Indicates that this product comes with Free Flow Cross Valve for fast drainage.	Curas created symbol
121.		Clear	Indicates that this product is clear (transparent) in colour.	Curas created symbol
122.		White	Indicates that this product is white in colour.	Curas created symbol
123.		Blue	Indicates that this product is blue in colour.	Curas created symbol
124.		Product capacity	Indicates that this product capacity is at volume indicated.	Curas created symbol
125.		Sugar free	Indicates that the product does not contain sugar.	Curas created symbol
126.		Lemon flavoring	Indicates that the product is with lemon flavor.	Curas created symbol
127.		Do not use on the eyes	Indicates that the product should not be used on the eyes.	Curas created symbol
128.		Approximate Graduation Mark	Indicates that this product has a graduation mark that showed approximate volume marking.	Curas created symbol
129.		Actual Graduation Mark	Indicates that this product has a graduation mark that showed actual volume marking.	Curas created symbol

SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
130.		Spout	Indicates that this product has lip projecting from a container, through which liquid can be poured.	Curas created symbol
131.		Maximum Width	Indicates that this product has a maximum width of 50mm.	Curas created symbol
132.		Heavy Metal Testing	Indicates that this product has been tested for heavy metal and the product pass the test according to the requirement.	Curas created symbol
133.		Maximum 5kg	Indicates that the load capacity is 5kg maximum.	Curas created symbol
134.		Tempered glass	Indicates that this product is made from tempered glass and has extra strength compared with normal glass.	Curas created symbol
135.		Auto zero	A function that forces the displayed frequency to return to <b>zero</b> when a pulse input is not received within a specified period of time.	Curas created symbol
136.		Load capacity	Indicates that the load capacity is less than or equal to the weight indicated.	Curas created symbol
137.		Approximate gross weight	Indicates the approximate gross weight of the product in kilograms.	Shipping & Packaging Mark Guide
138.		This way up	Indicates correct upright position of the transport package.	Shipping & Packaging Mark Guide
139.		Do not use hook	Indicates that hooks shall not be used for handling the transport package.	Shipping & Packaging Mark Guide

SYMBOLS GLOSSARY				
No.	Symbol	Symbol Title	Explanatory Text	Standard Reference
140.		Fragile, handle with care	Indicates that product can be broken or damaged if not handled carefully.	ISO 15223-1 Reference no. 5.3.1
141.		Do not stack	Indicates that vertical stacking is not allowed due to nature of the transport package.	Shipping & Packaging Mark Guide
142.		Recyclable	Indicates that the item can be recycled.	Shipping & Packaging Mark Guide
143.		Do not litter	Indicates that littering is prohibited.	Shipping & Packaging Mark Guide
144.		Recycle after use	Indicates that the product can be recycled after use.	Shipping & Packaging Mark Guide
145.		Suffocation hazard	Indicates that the product might suffocate children and must be kept away from children.	Shipping & Packaging Mark Guide
146.		Suffocation hazard	Indicates that the product might cause suffocation.	Shipping & Packaging Mark Guide
147.		Non-DEHP	Contains less than 0.1% Phthalates—DEHP	BS EN 15986:2011 EN 15986:2011(E) Annex B
148.		Type BF Part	To identify a type BF applied part. (Electrically connected to Patient but not directly to heart) Complying with IEC 60601	IEC 60417 Reference #5333 FDA Recognition #5-102
149.		EU Importer address	Indicates the medical device EU Importer.	ISO 15223-1