
































































Product Symbol Definition














Symbols	English (EN)
	European conformity
	CE marking directive, 93/42/EEC on medical devices
	Caution / Attention: See Instruction for Use (IFU)
	Consult instruction for use
	Do not re-use
	Contains or presence of natural rubber latex
	Latex-free
	Fragile, handle with care
	Keep dry
	Use by
	Date of manufacture
	Batch code
	Catalogue number















	Sterilized using ethylene oxide
	Sterilized using irradiation
	Non-sterile
	Do not use if package is damaged
	Manufacturer
	Authorized representative in the European Community
	Temperature limitation
	Keep away from sunlight
	Phthalate (DEHP)-containing device
	Phthalate (DEHP)-free
EN 455 (Parts 1, 2 & 3)	European standard for medical gloves
AQL 1.5	Acceptable quality levels of 1.5 for pinholes
	European standard: Protective gloves against chemicals and microorganisms. Comply with penetration test.





	<p>European standard: Protective gloves against chemicals and microorganisms. Microorganism resistant</p>
 <p>0000</p>	<p>European standard: Protective gloves against mechanical risks</p>
 <p>MIN. 4 HOURS</p>	<p>Fluid resistance test for a minimum of 4 hours (water at 35°C±3°C)</p>
<p>AS/NZS 4011 : 1997</p>	<p>Australian/New Zealand Standard: Single-use examination gloves specification</p>
<p>PPE Cat. 1</p>	<p>Personal Protective Equipment, European Directive 89/686/EEC, Category I: Protects from minor risks only</p>
	<p>Polyvinyl chloride-free</p>
	<p>Do not flush</p>
	<p>Period-after-opening</p>
	<p>Powder</p>
	<p>Powder-free</p>
	<p>Macerator friendly</p>
	<p>Pyrogenic-free</p>
	<p>For inflation of balloon catheter only</p>
	<p>Not for injection</p>

	Safe for contact with food
 Not For Oily Food	Not for oily food
	Non-return valve-free
	Kink-resistant
	Needle free sampling port connector
	Soap resistant maximum of 2 hours
	Polymer coated
	Fully textured
	Non-return valve
 Stainless Steel	Stainless steel 304
	Autoclavable
	Medium bag system
	Non-sterile inner bag system
	Single-patient-use system

	Perfume-free
	Colour free
	Low lint
	Lint-free
	High-density polyethylene
	Low-density polyethylene
	Polypropylene
	Slit cut
	Flexi-Funnel™
	Chlorine-free
	Antimicrobial
	Hypoallergenic adhesive
	Do not eat

	Keep out of reach from children
	Do not throw into basin
	Load capacity
	Light duty
	Medium duty
	Z folded for easy dispensing
	Safe for use in the freezer
	Star-shaped kink-resistant tube
	Safe for microwave
	Safe for oven
	Wash at 30°C
	Half folded for easy dispensing
	Formaldehyde-free

	High capacity cross valve
	Sampling port connector
	Open
 Z + 1/2 FOLD SINGLE DISPENSING	Z+1/2 folded for easy dispensing
	Chlorinated polyethylene
	Free flow cross valve
	White
	Clear
	Product capacity
	This way up
	Do not use hook
	Do not stack
	Recyclable
	Do not litter

	Recycle after use
	Suffocation hazard, keep away from children. This product is not a toy.
	Suffocation hazard
	Approximate gross weight

CUR/SYM/140715/00