

Weights Chart Guide

For Hygiene Pulp™ and C3 Emptying Bag w/SAP

This chart indicates the weight of Curas Hygiene Pulp™ products and C3 Emptying Bag in gramme. To weigh liquid in a Curas disposable Hygiene Pulp™ utensil, subtract from the total weight as follows:

REF	Product Description	Average Weight (g)	Tolerance +/- (g)
08000	Bedpan Liner, 2L, Large	70	5
08001	Slipper (Fracture) Pan Liner, 1.7L, Medium (Midi)	65	2
08002	Slipper (Fracture) Pan Liner, 2L, Large (Maxi)	74	2
08009	Bedpan Liner, 2L, Large, Non-Pressed	70	3
08026	2L Bowl w/Hygiene Control Bag	60	2
08027	3L Bowl w/Hygiene Control Bag	67	2
08028	4L Wash Bowl w/Hygiene Control Bag	107	5
08030	4L Wash Basin	73	5
08031	1L General Purpose Bowl	32	2
08032	3L General Purpose Bowl	67	2
08033	2L General Purpose Bowl	53	3
08034	Measuring Jug, 1L	40	2
08035	Sampling Cup, 200ml	11	2
08036	Ear Irrigation Cup (Noots Tank), 300ml	15	2
08040	Mid-Stream Urine Sampler, 1L	45	3
08046	Ear Irrigation Cup (Noots Tank), w/SAP by Curas™, 300ml	23	2
08050 / 08057	Urine Bottle, Round, 1L	50	5
08051 / 08055	Urine Bottle, Square, 850ml	45	5
08052	Urine Bottle, Female, 800ml	60	2
08056 / 08058	Urine Bottle, Round, 1L, w/SAP by Curas™ & Lid	62	5
08060	Urine Bottle, 1.5L, w/SAP by Curas™ & Lid, X-Large	92	5
08068/08070	Kidney Tray, 700ml	23	2
08071	Kidney Tray, 700ml BASIC	22	2
08073	1L General Purpose Bowl BASIC	22	2
08100	SAP by Curas™, Fluid Absorber, 6g Sachet	6	1
08120	Lid for Bedpan Liner	10	1
08121	Lid for Midi Slipper Pan Liner	13	1
08130	Lid for 1L General Purpose Bowl	7	1
08131	Lid for 2L General Purpose Bowl	13	1
08132	Lid for 3L General Purpose Bowl	12	1
08134	Lid for Jug	7	1
08135	Lid for Sampling Cup	2	1
08150	Lid for Round Urine Bottle	3	1

* All above weights exclude lid except REF 08056 & REF 08058.

Example:

Weight of urine bottle (REF 08050) with urine
*Weight of urine bottle (REF 08050) *see above chart*
Weight of urine
Volume of urine

252g
- 50g
<hr style="border: 0; border-top: 1px solid black;"/>
202g
≈202ml

For further info, kindly visit www.curas.com or contact us at info@curas.com.

Weights Chart Guide

For Hygiene Pulp™ and C3 Emptying Bag w/SAP

This chart indicates the weight of Curas Hygiene Pulp™ products and C3 Emptying Bag in gramme. To weigh liquid in a Curas disposable Hygiene Pulp™ utensil, subtract from the total weight as follows:

REF	Product Description	Average Weight (g)	Tolerance +/- (g)
08190	Graduated Stick for Male Urine Bottles, M, L and XL	5	1
08202	Deep Tray, S	16	2
08203	Deep Tray, M	70	3
08204	Deep Tray, L	30	2
08233	Procedure Tray w/PE, M	16	2
08568 / 08570	Kidney Tray, 700ml	24	2
08600	Urine Bottle, Round, 1L BASIC	30	2
08602	Bedpan Liner, 2L, Large BASIC	51	4
08603	Slipper (Fracture) Pan Liner, 1.7L, Medium (Midi) BASIC	44	4
09602	SAP by Curas™, Fluid Absorbent Pad, M size	12	1
09093	Bedpan Support, Heavy Duty	140	10
09094	Bedpan Support, Bariatric, Extra Heavy Duty	185	10
20460	C3 Emptying Bag w/SAP by Curas™, 2L	61	5
20465	C3 Emptying Bag w/SAP by Curas™, 1.5L	42	5
20467	C2 Short-term High Risk Drainage Bag w/SAP by Curas™, 2L	68	5
20469	C3 Emptying Bag w/SAP by Curas™, 2L, Non-DEHP	60	5
20470	C3 Emptying Bag w/SAP by Curas™, 2L, DEHP Free	41	5

* All above weights exclude lid except REF 08056 & REF 08058.

Example:

Weight of urine bottle (REF 08050) with urine
*Weight of urine bottle (REF 08050) *see above chart*
Weight of urine
Volume of urine

252g
- 50g
<u>202g</u>
≈202ml